

# Automatic Slide Stainer

## ASS 3700 Series



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

Automatic Slide Staining Machines not only provides high degree of flexibility and high specimen throughput in Histopathological Laboratories but also maintains maximum **Quality and Safety Standards**

**ASS-3700 Series AUTOMATIC SLIDE STAINING MACHINE (12 STATIONS):**

A simple, sturdy and economical model of automatic Staining Machine with 12 stations, recommended for hematology and also for staining of blood films and bone marrow spread etc.

►► **Key Features :**

- Compact Bench Top Model
- Perfect H&E Stainer for Routine Application
- High Specimen Throughput
- Minimum Reagent Exposure
- Easy to clean
- Simple to operate
- Highly recommended for Hematology and also for staining of blood Films and Bone marrow spreads etc.

►► **Technical Specifications :**

<b>Operation</b>	Automatic
<b>Staining Troughs</b>	12
<b>Troughs Capacity</b>	12 Slides
<b>Troughs Valume</b>	450 ml
<b>Intel Outlet Trough</b>	1
<b>Slide Carriers</b>	1
<b>1 Hr Calibrated Timing Discs</b>	6
<b>½ Hr Calibrated Timing Discs</b>	×
<b>Notch Cutting Plier</b>	√
<b>Dust Cover</b>	√
<b>Dimensions in mm</b>	560x220 V@50/60 Hz
<b>Power Supply</b>	110 / 220 V @ 50/60 Hz
<b>Timing Sequence</b>	12 Stage Programmble
<b>Timing Step</b>	1 Minute

▶▶ **HPLC Servicing, Validation, Trainings and Preventive Maintenance :**

- HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of HPLC’s Trainings.
- AMC’s/CMC : We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



▶▶ **About Analytical Technologies**

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

▶▶ **Our Products & Technologies**

						
UV/VIS Spectro 2080+ Double Beam	Infra FTIR	Optima Gas Chromatograph 3007	Optima Gas Chromatograph 2979 Plus	Flash Chromatograph	Atomic Absorption Spectrophotometer	Liquid Partical Counter
						
Optical Emission Spectrophotometer	DSC/TGA	Semi Auto Bio Chemistry Analyzer	HEMA 2062 Hematology Analyzer	Micro Plate Reader/Washer	URINOVA 2800 Urine Analyzer	Total Organic Carbon 3800
						
Fully Automated CLIA	NOVA-2100 Chemistry Analyzer	PCR/Gradient PCR/ RTPCR	TOC Analyzer	Laser Particle Size Analyzer	Ion Chromatograph	Water purification system

▶▶ Regulatory compliances



▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶ Reach us @



HPLC Solutions    MultipleLabs    Analytical Bio-Med    Analytical Distributors    Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.