



Blood and Infusion Warmer

Features:

- Delivers blood and intravenous fluid to the patient at Normothermic temperature at wide range of flow rates from gravity flow rates to 5,000 ml/hr
- Dedicated disposable triple lumen tubing that eliminates patient line cool down of infusate and delivers warm fluids to patients
- Built in safety feature that is over temperature test button and alarm test button

Technical Specification:

Physical		Dimensions
Height, Overall		24.1 cm (9.5 inches)
Width, Overall		21.0 cm (8.3 inches)
Depth, Overall		17.8 cm (7.0 inches)
Weight, Dry		3.5 Kg (7.6 lbs)
Weight, Wet (with recirculating solution)		5.0 Kg (11.0 lbs)
Weight, Shipping		3.6 Kg (7.95 lbs)
Recirculating Solution Capacity		1.4 L (0.37 gallons)
Maximum Height on I.V. Pole		107 cm (42 inches)
Environmental	Temperature	Humidity [%]
Operation	10°C to 45°C	10 to 95
Transportation	-18°C to 60°C	5 to 90
Storage	-18°C to 60°C	5 to 90
Thermal	Temperature	
Temperature Set Point	41.9oC ± 0.1oC	
Over Temperature Set Point	43.1oC	
Electrical	Type	
MAINS Power Input:		
230V	230VAC, 50/60 HZ, 1.5 Amps	
MAINS Auxiliary Supply Power		
Output:		
230V	230VAC, 50/60 HZ, 0.6 Amps	
Electrical	Type	
Protection Against Electrical Shock	Class 1 Equipment, Type BF	
Mode of Operation	Continuous	
Type of Current	Alternating	
Ingress Protection Rating	IPX1	
Performance		
Recirculating Solution		
Temperature	Recirculating solution temperature reaches 37oC from ambient in about 4 minutes	

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specification of the catalogue without prior notice