

USM-3008

Color Doppler Ultrasound Machine



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

▶▶ **High-end image within easy reach**

- Compound Imaging
- Speckle Reduction Imaging
- Tissue Harmonics Imaging
- New intelligent speckle noise suppression imaging technology



High light

High Resolution

High-LCD Monitor

High Sensitivity

Large Touch Screen

▶▶ **Main Features:**

- iBank database
- 4 probes connected
- Intelligent operation panel
- High quality real-time 4D technology
- 19 inch medical monitor + 10.4 inch touch monitor

Built-in 128G solid state drive, External USB storage

Rich clinical measurement software, apply to various department

Imaging Modes: 2D, 3D, 4D, Color/PW/CW/Power/Directional Color Power Doppler, Tissue Doppler, Free Steering M-Mode.

Probe Optional: Phased array probe, Convex probe, micro-convex probe, 4D volume probe, trans-vaginal probe, linear probe, trans-rectal probe.

▶▶ **Built-in integrated ultrasound workstation:**

- Direct patient file management
- Support image moving, static storage
- Input/Output: VGA, 4 USB Ports, HDMI, RJ-45, Built-in speakers
- Support for color and PW, three simultaneous display



Abdomen



Cardiology



Small Organ



Cavity Probe



Trans-Rectal Probe




Volume Probe

▶▶ **Super Good Images:**

High resolution cardiac software analysis with professional CW function

Application on: Cardiac, 4D Abdominal, OB/GYN, Urology, Vascular, Pediatric, MSK.



Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXIFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide