


COMFORT - DENT

DENTAL CHAIRS



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

COMFORT - TESLA



►► COMFORT- TESLA PROGRAMMABLE DENTAL CHAIR :

- Multiple Programmable Electrically • Operated Dental Chair, • (P1, P2, P3 Gargle and Synchroized programme), • Chair movement can be Operated by finger touch
- (capacitive touch) from instrument tray, and chair can be operated by feather touch from assistant side and also from multi foot control.

►► Dental Unit consisting of :

- LED Light head with SENSOR operating • Glass spittoon/tumbler with auto water system
- Adjustable, up/down instrument tray with s.s.tray • Attached Round Glass tray
- X-Ray viewer • 2.Nos Three-way syringe • 2.Nos Air-roter connection • Micro-motor connection
- Suction - One medium motorized suction and one Air Suction
- Disinfection Bioflush System

COMFORT - JET 3200



►► COMFORT- JET 3200 BASE-LESS DENTAL CHAIR MOUNT UNIT :

- Electrically Operated Chair mount unit with zero position
 - Imported CE Marked Motor for chair movement
 - PU moulded arm rest
 - Wide Instrument tray
 - Glass /Ceramic spittoon & cup filler with auto timer system
 - Imported halogen light with 2 intensity
 - X-ray film viewer
 - 2 No. 3 way syringe
 - 2No. Air-rotor point
 - Micro-motor point
 - Suction - One medium motorized suction and one Air Suction
- Optaonal : Scalar, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - ELENTRA



►► COMFORT- ELENTRA DENTAL CHAIR MOUNT UNIT :

- Electrically Operated Chair mount unit with zero position
 - Imported CE Marked Motor for chair movement
 - PU moulded arm rest
 - Wide Instrument tray
 - Glass /Ceramic spittoon & cup filler with auto timer system
 - Imported halogen light with 2 intensity
 - X-ray film viewer
 - 2 No. 3 way syringe
 - 2No. Air-rotor point
 - Micro-motor point
 - Suction - One medium motorized suction and one Air Suction
- Optaonal : Scalar, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - JET 3100



Woodpecker Scaler



Micro Motor

►► COMFORT- JET 3100 BASE-LESS DENTAL CHAIR MOUNT UNIT :

- Electrically Operated Chair mount unit with zero position
 - Imported CE Marked Motor for chair movement
 - PU moulded arm rest
 - Wide Instrument tray
 - Glass /Ceramic spittoon & cup filler with auto timer system
 - Imported halogen light with 2 intensity
 - X-ray film viewer
 - 2 No. 3 way syringe
 - 2No. Air-rotor point
 - Micro-motor point
 - Suction - One medium motorized suction and one Air Suction
- Optaonal : Scaler, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - PEAK



LED Light Cure



Oil Free compressor

►► COMFORT-PEAK SUPER DELUX DENTAL CHAIR MOUNT UNIT :

- Electrically Operated Chair mount unit with zero position
 - Imported CE Marked Motor for chair movement
 - PU moulded arm rest
 - Wide Instrument tray
 - Glass /Ceramic spittoon & cup filler with auto timer system
 - Imported halogen light with 2 intensity
 - X-ray film viewer
 - 2 No. 3 way syringe
 - 2No. Air-rotor point
 - Micro-motor point
 - Suction - One medium motorized suction and one Air Suction
- Optaonal : Scalar, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - DELUX



Auto Clave



Dental X-ray machine

►► COMFORT- DELUX ELECTRICALLY OPERATED DENTAL CHAIR MOUNT :

- Dental Unit consisting of • Halogen light • X-Ray viewer • Three Way Syringe with straight tubing
- Air Roter connection with straight tubing • Micro motor connection
- Pneumatic High & Low suction.

Optional: Scalar, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - VISION



Straight Handpiece



Contra Angle Handpiece



Air-Rotor Handpiece



Doctor Stool

►► COMFORT- VISION FOLDABLE DENTAL CHAIR MOUNT UNIT :

- Electrically Operated Chair mount unit with zero position • Imported CE Marked Motor for chair movement • PU moulded arm rest • Wide Instrument tray • Glass /Ceramic spittoon & cup filler with auto timer system. Imported halogen light with 2 intensity • X-ray film viewer
- 2 No. 3 way syringe • 2No. Air-rotor point • Micro-motor point
- Suction - One medium motorized suction and one Air Suction
- Optional : Scalar, Light Cure, Micro Motor with Hand piece, Air-Rotor Hand piece.

COMFORT - KID PEDO



►► COMFORT- KID PEDO DENTAL CHAIR :

- Electrically Operated Dental Chair, • Chair movement can be by feather touch control fitted with instrument tray, assistant side and also by multi foot control,

Dental Unit consisting of:

- Halogen Light with two intensity. Glass spittoon/tumbler with auto water system
- Adjustable, up/down instrument tray with s.s.tray . Attached Glass round tray. X-Ray viewer
- 2.Nos Three-way syringe . 2.Nos Air Roter connection • Micro-motor connection
- Suction - One medium motorized suction and one Air Suction

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide