

Digital Overhead Stirrer 3000



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Specifications

Dimension	420*270*200 mm
Packing size	48*33*25 cm
Weight	12 kg
Max. stirring quantity	20L (H ₂ O)
Motor rating input,	40 W
Motor rating output	50 W
Speed range	40-2500 rpm
Speed display accuracy	±3 rpm
Speed display	LCD
Speed display resolution	1 rpm
Max. torque	40 Ncm
Torque display	LCD
Overload protection	LCD light flash, auto stop
Motor protect	LCD light flash, auto stop
Viscosity max. mPas	10000
Chuck range diameter	0.5-13 mm
Rod size dia.	12x230 xl, mm
Vlotage	100-240 VAC
Frequency	50/60 Hz
Power	50 W
Protection class acc	DIN/EN 60529 IP42
Temperature	5-40 °C
Permission relative humidity	80 %
Remote control	RS232 interface

►► Servicing, Validation, Trainings and Preventive Maintenance :

Servicing	: We have team of service engineers who can attend to any make of instrument promptly @the most affordable cost.
Trainings	:We also take up preventive maintenance to reduce downtime of instrument's Trainings.
AMC's/CMC	: We offer user training both in-House and at customer sites on instrument principles, operations, troubleshooting.
Validations	: We have protocols for carrying out periodic Validations as per GLP/ GMP/USFDA norms.
Instruments	: We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



►► About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

►► Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

►► Regulatory compliances



►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

►► Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.