

## DIGITAL VISCOMETER



The Digital DV-E is the low cost solution for those customers wanting a calculating, digital, basic viscometer. This is a good upgrade choice for those previously using dial reading instrument.

### ►►► Specification:

- Low cost digital viscometer
- Easy to use – Flip a switch  
Turn a knob
- 18 Speeds for complete range capability
- No calculations required – Direct reading of viscosity in cP or mP
- Complete with appropriate spindles, viscometer stand, guard leg and carrying case.

|                               |                  |
|-------------------------------|------------------|
| Repeatability:                | :± 0.2%          |
| Accuracy:                     | :± 1.0% of range |
| Range Minimum                 | :1+ Range        |
| Maximum                       | :2M.             |
| Speed                         | :3.100 RPM       |
| In Speed Number of increments | :18              |
| Spindles                      | :Included.       |



## ▶▶▶ Regulatory compliances



## ▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## ▶▶▶ Reach us @



HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide