

## Incubator Bacteriological



Standard double wall fabrication. Inner chamber made of richly anodized aluminum or highly polished stainless steel sheet and outer made out of thick mild steel finished in white stoving enamel/powder coated paint with mat finished colour combinations. Double wall door with double glass window for observation in the chamber. Temperature range 5°C above ambient to 80°C  $\pm$  1°C controlled by a thermostat. The equipment is workable on 220V AC 50 Hz single phase.

### (Chambers size in mm/inches)

(H x W x D) Approx.	Shelves	No. of Capacity (Approx)
1. 300 x 300 x 300	(12" x 12" x 12")	2 28 ltrs.
2. 350 x 350 x 350	(14" x 14" x 14")	2 45 ltrs.
3. 400 x 400 x 400	(16" x 16" x 16")	2 65 ltrs.
4. 450 x 450 x 450	(18" x 18" x 18")	2 95 ltrs.
5. 600 x 450 x 450	(24" x 18" x 18")	2 125 ltrs.
6. 600 x 600 x 600	(24" x 24" x 24")	2 224 ltrs.
7. 900 x 600 x 450	(36" x 24" x 18")	2 252 ltrs.
8. 900 x 600 x 600	(36" x 24" x 24")	2 336 ltrs.
9. 900 x 900 x 900	(36" x 36" x 36")	3 480 ltrs.
10. 1200 x 600 x 840	(48" x 24" x 34")	3 600 ltrs.

### Optional:

- 1) Air circulating fan.
- 2) Timer 0-24 hours
- 3) Imported Thermostat "EGO" / "JUMO" German with the accuracy of  $\pm$  0.5°C.
- 4) Digital temp indicator-cum-Controller.
- 5) Microprocessor PID digital temp. Indicator-cum-controller.

## Regulatory compliances



## Corporate Social Responsibility



Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.



HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide