

Ice flaker

Features:

- Fully Automatic microprocessor control
- Stainless Steel construction
- Continuous Ice Flakes Output
- CFC free compressor
- Low water level Indication
- Overload protection
- Noiseless operation
- Compact Design and low maintenance



Technical Specification:

Capacity	300kg/day
Storage Bin	ABS Plastic
Cooling Method	Air
Working temperature	Below 30 °C for optimum working
Construction	Stainless Steel
Noise Level	<52 dB
Water Supply	3/4 inch BSP
Input Power	180W

Regulatory compliances



Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for an award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFIFY human minds, souls and bodies by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

Reach us @



HPLC Solutions	MultipleLabs	Analytical Bio-Med	Analytical Distributors	Analytical Foundation (Trust)
Corporate & Regd. Office: Analytical House, # E67 & E68, Ravi Park, Vasna Road, Baroda, Gujarat 390 015. INDIA	T: +91 265 2253620 +91 265 2252839 +91 265 2252370 F: +91 265 2254395	E: info@hplctechnologies.com info@multiplelabs.com info@analyticalgroup.net info@analyticalbiomed.com	W: www.ais-india.com www.analyticalgroup.net www.hplctechnologies.com www.multiplelabs.com	Sales & Support Offices: across the country : Distributors & Channel partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.