

# NOVA-2600

Auto Dry Chemistry Analyser



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## ►► Specifications

|                       |   |
|-----------------------|---|
| Sample type:          | Whole blood, serum, plasma  |
| Sample volume:        | 100ul   |
| Test time :           | About 10 minutes  |
| Parameters:           | Up to 17 reportable parameters in a test, 32 parameters configured into 9 profiles.   |
| Test Profile:         | ALB, ALKP, ALT, AMYL, AST, BUN, Ca, CHOL, CK, CREA, GGT, GLU, LAC, LDH, LIPA, NH3, PHOS, TBIL, TP, TRIG, URIC, GLOB, UCRE, UPRO, UPC Ratio, Mg, OSMOL |
| Reagent disc:         | Disposable, self-contained reagent in dry form  |
| Reagent information:  | Built-in QR code  |
| Interface:            | WLAN, One USB ports, 1 Ethernet interface, support LIS  |
| Display:              | 3.5 inch color LCD touch screen   |
| Temperature Accuracy: | 37°C ± 0.1°C  |
| Working temperature:  | 10-30°C Humidity: 40-85%  |
| Power:                | AC100V-240V, 50-60Hz, DC 24V  |
| Printout :            | Built-in thermal printer(optional), chemistry management platform   |
| Storage capacity:     | More than 50000   |
| Weight :              | About 2.5KG   |
| Dimension:            | 175 (H) * 125 (W) * 210(L) MM   |

## ►► Features:

- Smaller--only 2.5kg
- Accurate--calibration, quality control
- Faster--up to 17 results in just 7-13min
- Easier--3 steps operation



## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies




UV/VIS Spectro 2080+ Double Beam



Infra FTIR



Optima Gas Chromatograph 3007



Optima Gas Chromatograph 2979 Plus



Flash Chromatograph



Atomic Absorption Spectrophotometer



Liquid Partical Counter



Optical Emission Spectrophotometer



DSC/TGA



Semi Auto Bio Chemistry Analyzer



HEMA 2062 Hematology Analyzer



Micro Plate Reader/Washer



URINOVA 2800 Urine Analyzer



Total Organic Carbon 3800



Fully Automated CLIA



NOVA-2100 Chemistry Analyzer



PCR/Gradient PCR/ RTPCR



TOC Analyzer



Laser Particle Size Analyzer



Ion Chromatograph



Water purification system

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide