

Trinocular Microscope



Specification:

- ▶ Frame: All aluminum pressure Diecasted body.
- ▶ Paint: P.U. coating in light cream shade.
- ▶ Nose piece: Quadruple nosepiece.
- ▶ Binocular head: 360° rotate able and 45° inclined head with all anti-fungus coated prisms.
- ▶ Focusing system: Coarse and fine motion moving on rack and pinion
- ▶ Mechanical stage: Cxl mechanical stage for smooth functioning
- ▶ Condenser: Condenser Bright field sub stage condenser having 1.2 N.A.
- ▶ Objectives: 4x, 10x, 40x, and 100x.
- ▶ Eye piece: WF-10X (Pair)
- ▶ Illumination: 6V, 20w halogen lamp needs 220v (50Hz), 110(60Hz) optional supply.

Analytical Technologies, USB Microscope camera:

- ▶ USB Camera: 5MP USB high Resolution Color cooled camera with Photo

Lens, Microscope adapter specially recommended for fluorescence Microscope

- ▶ Sensor type: cooled CCD
- ▶ Sensor format: 2/3", with software

Regulatory compliances



Corporate Social Responsibility



Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Workshops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFIFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide