

Urine Analyzer-UA-2100



ENDLESS CHALLENGE
ENDLESS IMPROVEMENT

EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

ATL Urine Analyzer-UA-2100



Test Items:	Glucose. Bilirubin. Ketone. Specific Gravity. Blood.pH. Protein. Urobilinogen Nitrite. Leukocyte. Vitamin C
Test Method:	Dry strip color-developing method
Test Principle:	Reflective photometry with dual-wavelength
Light Source:	Cold lamp with high brightness and long lifetime
Wavelength:	520nm,640nm,660nm
Test Mode:	Single and continuous test mode are selectable
Throughput:	120 Strips/Hour
Test Slot:	Using Sliding Slot, Easy to clean and replace, Lower the Carryover
Test Strip:	Provide original strips
Display:	LCD with high resolution
Language:	Chinese and English are selectable
Memory:	Saving more than 2000 samples results, and auto recover the results when power off during working
Interface:	RS-232 232 serial porte can be connected to external printer and computer
Printer:	Inner thermal pñnter can be connected to external printer
Report Format:	Half quantitative and qualitative
Result Unit:	SI unit and routine unit are selectable
Status Check:	Auto checking the status during power on and off auto indentify the alarm.
Working Environment:	Temperature: 10 C-30U ;Hurnidity: 20%RI-1-80%RH
Voltage:	85 264V,Auto adaption according to the input voltage
Power:	<80W

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/ GMP/USFDA norms.
- Instruments** :Instruments :We offer instruments / Renting Services Modules like pumps, detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies

						
UV/VIS Spectro 2060+ Double Beam	FTIR	Gas Chromatograph 3000	Gas Chromatograph 2979 Plus	Flash Chromatograph	Atomic Absorption Spectrophotometer	Liquid Partical Counter
						
Optical Emission Spectrophotometer	DSC/TGA	NOVA 2020 plus Automated Bio Chemistry Analyzer	HEMA 2020 Hematology Analyzer	Micro Plate Reader/Washer	URINOVA 2800 Urine Analyzer	Total Organic Carbon
						
Fully Automated CLIA	NOVA Basic Semi-Auto Chemistry Analyzer	PCR/Gradient PCR/ RTPCR	Blood Gas Analyzer	Random access Analyzer for immunoassay Proteins & clinical chemistry	Semen Analyzer	Water purification system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.