



Digital Water Baths with 20L DWB-3001 Series



Made available in options of digital model and thermostatic model options

Suitable for holding 2, 4, 6, 8, racks

Available with thermostatic based temperature control support

Stirring option available with/without stirring arrangement with support of FHP stirrer (optional)

Comes with glass mercury thermometer based temperature display provision

Featuring mineral wool based insulation support

Perforated cover of stainless steel based diffuser shelf finish

Specification

Capacity	20L
Heating Load	4.0KW
Temperature Controller	Thermostatic
Stirring (Optional)	Available With or without stirring Arrangement with the help of FHP stirrer
Temperature Range	Ambient +5°C - 97°C +0.5°C
Temperature Display	Glass Mercury Thermometer
Insulation	Mineral Wool
Diffuser Shelf	Perforated Cover of stainless steel
Electrical Supply	220/230V AC, 50/60Hz

*Note: Company Can Custom/Tailor make as per the customer orders/requirements.

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for an award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and bodies by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

▶▶▶ Reach us @



Analytical[®]
Technologies Limited

HPLC Solutions | MultipleLabs | Analytical Bio-Med | Analytical Distributors | Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide