



HOT AIR OVEN HAO-3001 Series



Technical Specification

Inner Chamber Size	450x450x450 mm
Volume	125 litre
Temperature variation (Time)	+/-0.5°C
Temperature deviation	+/-0.5°C
Readability/ Set ability	0.5°C
Temperature range	25-300°C
Degree	above ambient to 250 degree C
Sensor thermocouple	Type J
Controller	Nano technology based ergonomic controller
Display	220 YAC, 1 Phase, 50Hz
Watts	1750 - 2000
Trays	2 or 3

*Note: Company Can Custom/Tailor make as per the customer orders/requirements.

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) founded for the purpose of:

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various fields of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Workshops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFIFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Comps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplc technologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide